



Welcome to the **New World** of **Immediate**

PALTOP System
 Four-Tier Advantage



Paltop designed a 3rd generation implant taking the benefits found in both tapered implants and cylindrical implants. Furthermore, its design enables the same drills to be used regardless of implant length. This marriage creates easy and safe insertion along with bone condensation for stability; providing clinicians the freedom to treat patients with an immediately or delayed application, all using the same system.

Implants have been developed "from biology to design" enabling:

- 1 Platform shifting -to minimize bone resorption
- 2 Micro threads to reduce bone stresses and minimize bone resorption
- 3 Cylindrical profile for bone compression and initial stability for immediate loading
- 4 Cutting threads and tapered design enable smooth and easy insertion
- 5 Passive apex for safe insertion in challenging anatomical positions

These benefits translate into providing clinicians with greater safety and predictability of results while maintaining flexibility to serve each patient's unique set of needs not compromising costs and simplicity.

3rd Generation Implant Design

