


1.25 Drivers

NEW ARRIVALS



60-70003

**Hex Driver for
Ratchet Ø1.25mm,
Long**

This driver connects to the ratchet without the need for an adaptor.

60-70164

**Hex Key Ø1.25mm
L6.5mm**

This driver's extra short length assists with hard to reach areas. It is used directly in the contra angle.

60-70165

**Hex Key Ø1.25mm
L23mm**

This driver's extra long length assists with tall restorations. It is used directly in the contra angle.

60-70005

**Manual Hex Driver
Ø1.25mm, Short**

This driver's extra short length assists with hard to reach areas. Now with a new design for a better grip.