



Gold Base Castable Abutments

Superior option for fabricating implant-level custom cast restorations.

Paltop Gold Base Castable Abutment may be used for casting individual abutments and screw-retained crowns/bridges directly to the implant connection.

These abutment assemblies consist of a gold alloy interface implant connection for high accuracy and excellent sealing, and a burn-out plastic cylinder for casting.

Paltop offers Gold Base Castable Abutments for Advanced Standard Platform, Advanced Narrow Platform and Conical connections. For each Platform, an Engage and a Non-Engage connection is available for multiple implant bridges.



**Standard
Platform
Engage**
46-70020



**Standard
Platform
Non-Engage**
46-70022



**Narrow
Platform
Engage**
46-70019



**Narrow
Platform
Non-Engage**
46-70021



**Conical
Engage**
46-72001



**Conical
Non-Engage**
46-72002