

NEW & IMPROVED! **CHROME COBALTBASE**

New option for fabricating implant-level custom cast restorations.

Paltop Chrome Cobalt Base Castable Abutment may be used for casting individual abutments and screw-retained crowns/bridges directly to the implant connection.

These abutment assemblies consist of a chrome cobalt interface implant connection for high accuracy and excellent sealing, and a burn-out plastic cylinder for casting.

Paltop offers chrome cobalt Base Castable Abutments for Standard Platform, Narrow Platform, Wide Platform and Conical connections. For each Platform, an Engage and non-engage platform is available to manage single tooth and bridge application.

BENEFITS OF OUR IMPROVED DESIGN:

- **Now offering 4 platforms: NP, SP, WP & PCA**
- **New narrow base to allow greater flexibility in developing the emergence profile.**
- **Tighter and stronger connection between castable sleeve and chrome cobalt base.**



**Standard
Platform
Engage**

40-70077



**Standard
Platform
Non-Engage**

40-70177



**Narrow
Platform
Engage**

40-70078



**Narrow
Platform
Non-Engage**

40-70178



**Wide
Platform
Engage**

40-70174



**Wide
Platform
Non-Engage**

40-70175



**Conical
Engage**

40-72021



**Conical
Non-Engage**

40-72022