

PALTOP

**NEW
ARRIVALS**


SHORT STEP-DRILL KIT





cat #- 60-70150

Paltop's short step-drills, now provided in a new small kit.

The Short Step Drill Kit is useful when there is limited interarch space for drilling, it has 4 stepped twist drills:

1.5/2.0, 2.0/2.4, 2.4/3.2, 3.2/3.8



1.5/2.0 mm Short Drill		60-70047
2.0/2.4 mm Short Drill		60-70049
2.4/3.2 mm Short Drill		60-70051
3.2/3.8 mm Short Drill		60-70053