

# LOCATORS ABUTMENTS




For over 30 years, ZEST Anchors has been a global leader in the manufacturing of dental implant attachments.

The LOCATOR implant Attachment is the premier system for implant retained overdentures, with the benefit of combating the damage done by the improper seating of overdentures.

**Now available for PALTOP implants  
for the Advanced and the Conical lines.**

## LOCATOR® MALES

The unique Dual Retention innovation provides the LOCATOR Attachment with a greater retention surface area than ever before available with other attachments.



1.5 lbs



3 lbs



5 lbs

## LOWEST VERTICAL HEIGHT

The LOCATOR Male remains in static contact with the attachment while the Denture Cap, which is processed into the overdenture, has a full range of rotational movement over the male for a genuine resilient connection of the prosthesis without any resulting loss of retention.

Locator	ERA	EDS
3.17mm	4.85mm	6.22mm



Dual retention, pivoting action provides resiliency to maximize stability and longevity.



# LOCATORS

## FOR PALTOP IMPLANTS

<b>40-70082</b>	Locator for Paltop - Advanced - standard platform 1mm Cuff Height
<b>40-70083</b>	Locator for Paltop - Advanced - standard platform 2mm Cuff Height
<b>40-70084</b>	Locator for Paltop - Advanced - standard platform 3mm Cuff Height
<b>40-70085</b>	Locator for Paltop - Advanced - standard platform 4mm Cuff Height
<b>40-70086</b>	Locator for Paltop - Advanced - standard platform 5mm Cuff Height
<b>40-70087</b>	Locator for Paltop - Advanced - standard platform 6mm Cuff Height

<b>40-70088</b>	Locator for Paltop - Advanced - Narrow Platform 1mm Cuff Height
<b>40-70089</b>	Locator for Paltop - Advanced - Narrow Platform 2mm Cuff Height
<b>40-70090</b>	Locator for Paltop - Advanced - Narrow Platform 3mm Cuff Height
<b>40-70091</b>	Locator for Paltop - Advanced - Narrow Platform 4mm Cuff Height
<b>40-70092</b>	Locator for Paltop - Advanced - Narrow Platform 5mm Cuff Height
<b>40-70093</b>	Locator for Paltop - Advanced - Narrow Platform 6mm Cuff Height

<b>40-72011</b>	Locator for Paltop - Conical - 1mm Cuff Height
<b>40-72012</b>	Locator for Paltop - Conical - 2mm Cuff Height
<b>40-72013</b>	Locator for Paltop - Conical - 3mm Cuff Height
<b>40-72014</b>	Locator for Paltop - Conical - 4mm Cuff Height
<b>40-72015</b>	Locator for Paltop - Conical - 5mm Cuff Height
<b>40-72016</b>	Locator for Paltop - Conical - 6mm Cuff Height

<b>40-70081</b>	Locator Male Processing Package
-----------------	---------------------------------

