



WIDE PLATFORM 6MM WIDE IMPLANTS



The PALTOP Wide Platform 6.0mm implant has an outer diameter of 6.0mm and a 4.5mm hexagonal connection. It is designed for restoration of the posterior regions and molar extraction sites. The PALTOP Wide Platform has a matching wide prosthetic line.



ADVANCED CLASSIC

6.0	
6.0 mm	20-70022
8.0 mm	20-70023
10.0 mm	20-70024
11.5 mm	20-70025
13.0 mm	20-70026
16.0 mm	20-70027



ADVANCED PLUS

6.0	
6.0 mm	20-70022P
8.0 mm	20-70023P
10.0 mm	20-70024P
11.5 mm	20-70025P
13.0 mm	20-70026P
16.0 mm	20-70027P







DYNAMIC

6.0	
6.0 mm	21-70022
8.0 mm	21-70023
10.0 mm	21-70024
11.5 mm	21-70025
13.0 mm	21-70026
16.0 mm	21-70027

6MM WIDE IMPLANT DRILLING KIT

cat # 60-70114



60-70098	6.0 mm implant long final drill	
60-70121	3.8/4.3 Long step twist drill	
60-70123	4.3/5.3 Long step twist drill	
60-70143	Countersink for 5mm&6mm Implant	

WIDE PLATFORM PROSTHETICS


CONCAVE HEALING CAPS




STRAIGHT HEALING CAPS




FULLY CASTABLE ABUTMENTS




TEMPORARY TITANIUM ABUTMENTS




MULTI-UNIT



IMPRESSION COPINGS



TITANIUM BASES



STRAIGHT TITANIUM ABUTMENTS

